

Vitamin C Content of Foods

Experiment #11

Pre-lab Exercise

1. Vitamin C is considered to be unstable during cooking. What kind of chemical reaction do you expect vitamin C to undergo in solution (hydrolysis, oxidation, reduction)?
2. Show the chemical reaction you would expect for loss of ascorbic acid when heated in air (O_2).
3. The formula shown above (eq 1) can be used to calculate the amount of ascorbic acid (mg) in test solutions. If 12.5 ml of iodine solution is needed to titrate 10.0 ml of 1.0 mg/ml ascorbic acid standard solution and 5.0 ml of iodine solution is needed to titrate 10.0 ml of an ascorbic acid test solution. What is the concentration of ascorbic acid in the test solution?

4. In Experiment #7, The Chemistry of Carbohydrates, you tested for starch with iodine solution. What color change is observed when iodine solution is added to starch solution?

5. What is the recommended dietary allowance (RDA) for vitamin C for adults?

6. What disease results if a person becomes deficient in vitamin C? What are the symptoms of this disease?

7. Given the symptoms of vitamin C deficiency, what is one biochemical function for vitamin C in the body? You should describe the biochemical or metabolic function in which vitamin C is involved (see textbook).

Name _____

Section _____

Vitamin C Content of Foods

Experiment #11

Data & Report Sheet

Part B. Titration of Ascorbic Acid Standard Solution (1 mg/ml):

	Flask 1	Flask 2
Final level of I ₂ soln. in buret: (Read from top of buret down)	_____	_____
Initial level of I ₂ soln. in buret: (Zero is at top of the buret)	- _____	- _____
Total vol. of I ₂ soln. used:	_____	_____
Average volume used:	_____ ml	

B-1. When you mix 10.0 ml of ascorbic acid standard solution in the flask with 25 ml water, 2 ml of 6M acetic acid solution and 1 ml of starch solution, is it necessary to know the exact final volume of this solution to get an accurate determination of ascorbic acid in the titrations of test solutions? Explain. Hint: Do the volumes of water, acetic acid and starch solutions enter into the calculations to determine the amount of ascorbic acid? Does the volume of ascorbic acid solution added before titration influence how much iodine solution will be needed?

Part C. Determination of Ascorbic Acid in Test Solutions:

Titration Results	deionized H ₂ O		alkaline solution	
	1	2	1	2
Fin vol, ml				
Init vol, ml				
Tot used, ml				
Vit C, mg/ml				
% Original Vit C Remaining in Sample				

C-1. Give a brief explanation for any differences in the amount of ascorbic acid remaining in the heated test solutions, *i.e.*, what may cause differences between deionized water and alkaline conditions?

C-2. What effect does heating have on any disappearance of vitamin C?

Name _____

Section _____

Part D. Vitamin C in Beverages.

Beverage Samples: A) _____ mL of _____
Type of Beverage

B) _____ mL of _____
Type of Beverage

C) 5 mL of vitamin C tablet dissolved in 100 mL water.

Titration Results

	A	B	C
Fin Vol, ml			
Init vol, ml			
Total vol, ml			
Conc Vit C, mg/ml			
Vit C/serving (mg/ml x ml/serv)			
% RDA for Vit C (RDA = 60 mg)			

D-1. Check the labels for the juice or beverage containers and record whether vitamin C is listed for that beverage and indicate how much is supposed to be in the beverage if that information is given. If it's not given, do you think it should be?

D-2. Would you consider either of the beverage samples a "good" source of vitamin C in the diet? Give your own clarification of what a "good" source would be.

[Answer question on back of this page].

D-3. If you measured vitamin C in the solution prepared from a vitamin C tablet. Calculate how much vitamin C was in the tablet (show your work) and compare that with the amount claimed to be in the tablet given on the label. What may account for any significant difference between what you measured and what is claimed?

Give a brief summary of your own conclusions about the chemical stability of vitamin C and how this might affect the amount of vitamin C that may be found in foods, such as cooked vegetables.