

# The Chemistry of Carbohydrates


## Experiment #7

**Objective:** To determine the carbohydrate class of an unknown by carrying out a series of chemical reactions with the unknown and known compounds in each class of carbohydrates.


### Introduction

Simple sugars, starches and cellulose are organic compounds that have the approximate formula  $C(H_2O)_n$ , which accounts for the name carbohydrate (or hydrate of carbon) that is usually applied to this group of compounds. They are not truly hydrates of carbon but are polyhydroxy (alcohol) compounds that contain an aldehyde or ketone functional group. These functional groups give the carbohydrates some of their chemical properties that will be studied in this lab.

Simple sugars are called monosaccharides (one sugar), or disaccharides (2 sugars). Some monosaccharides are glucose, fructose, galactose, and xylose. Note that xylose is a pentose and fructose is a ketose.




Two common disaccharides are sucrose (table sugar) and lactose (milk sugar); sucrose is a combination of glucose and fructose linked together by their anomeric carbons to produce a nonreducing sugar (it does not reduce  $\text{Cu}^{2+}$ ), whereas lactose is a combination of galactose and glucose linked together by a  $\beta$ -1,4-glycosidic bond to produce a reducing disaccharide.



When many sugar molecules are linked together into a polymer, the resulting compound is called a polysaccharide. Starches and celluloses are polysaccharides. Amylose is a linear chain polymer of glucose, whereas amylopectin (a plant starch) and glycogen (an animal starch)

are branched polymers of glucose (see the text book for structures of the branched starches).




### Qualitative Tests for Carbohydrates

Reducing sugars are usually detected with ***Benedict's reagent***, which contains  $\text{Cu}^{2+}$  ions in alkaline solution with sodium citrate added to keep the cupric ions in solution. The alkaline conditions of this test causes isomeric transformation of ketoses to aldoses, resulting in all monosaccharides and most disaccharides reducing the blue  $\text{Cu}^{2+}$  ion to cuprous oxide ( $\text{Cu}_2\text{O}$ ), a brick red-orange precipitate. This solution has been used in clinical laboratories for testing urine.

***Barfoed's solution*** contains cupric ions in an acidic medium. The milder condition allows oxidation of *monosaccharides* but does not oxidize *disaccharides*. If the time of heating is carefully controlled, disaccharides do not react while reducing monosaccharides give the positive result (red  $\text{Cu}_2\text{O}$  precipitate). Ketoses do not isomerize with this reagent.

Carbohydrates are dehydrated in the presence of nonoxidizing acids to form furfural and hydroxymethylfurfural.



***Seliwanoff's reagent*** contains resorcinol in 6 M hydrochloric acid. Hexoses undergo dehydration when heated in this reagent to form hydroxymethylfurfural, that condenses with resorcinol to give a red product. Keto-hexoses (such as fructose) and disaccharides containing a keto-hexose (such as sucrose) form a cherry-red condensation product. Other sugars may produce yellow to faint pink colors.

***Bial's reagent*** contains orcinol (5-methylresorcinol) in concentrated HCl with a small amount of  $\text{FeCl}_3$  catalyst. Pentoses are converted to furfural by this reagent, which form a blue-green color with orcinol. This test is used to distinguish pentoses from hexoses.

**Iodine** forms a deep blue color in the presence of starch. Potassium iodide is added to the reagent solution in order to make the iodine more soluble in water. Some forms of starch may yield a greenish color. Simple carbohydrates (mono- and disaccharides) and cellulose do not cause any change in the orange-brown color of the iodine reagent.

### **Reagents and Materials**

Packet of carbohydrate unknown, 1% solutions of carbohydrate standards (glucose, fructose, xylose, sucrose, lactose, starch), Benedict's reagent, Barfoed's reagent, Selinwanoff's reagent, Bial's reagent, iodine in potassium iodide solution, testtubes, beaker for hot water bath, bunsen burner.

**\*\*\* Begin by setting up tubes for the enzymatic digestion of starch in Part F, steps 1-5.**

### **Preliminary Procedures - Do These Before Starting Part A.**

- Obtain an unknown carbohydrate and prepare a 1% solution of your unknown by dissolving 0.25 g in 25 mL of deionized water. Test your unknown carbohydrate solution in parallel with the known carbohydrate solutions in each of the following tests. If you work in pairs, each student in the group must have his or her own unknown.
- Set up a hot water bath using a 400 mL beaker that is about half filled with water and heat to near boiling. Do not fill more than half full.
- **Note:** Begin by setting up the mixtures for part F before proceeding further in part A.
- The purpose of this experiment is to demonstrate that different classes of carbohydrates can be distinguished from one another by specific chemical tests.

### **Part A: Test for reducing sugars with Benedict's reagent.**

Be sure to start Part F before proceeding.

1. Label 8 test tubes for each of the 7 test carbohydrates and water blank. Test carbohydrates are: 1-glucose; 2-fructose; 3-xylose; 4-sucrose; 5-lactose; 6-starch; and 7-unknown. Tube number 8 is for a water blank. This scheme is used for parts A, B and C.
2. Add about 3 mL of 1% carbohydrate test solution (those listed above) to each of the labeled test tubes. Add 3 mL of deionized water to tube #8 as a blank. If you are working in pairs, you will need an additional tube for the second unknown.

**Caution: Benedict's reagent is caustic, rinse thoroughly with water if you get this solution on your skin or clothing.**

3. Add 3 drops of Benedict's solution to each tube and mix thoroughly.

4. Place the test tubes in the hot water bath and note the color changes within 2 min as they are heated (the water bath should be hot, but not boiling, heat the tubes no more than 3 min).
5. Record observations on your report sheet.
6. What conclusions can you make about your unknown?

These solutions can be discarded in the sink.

**Part B: Test for mono- vs disaccharides with Barfoed's reagent.**

1. Empty, wash and rinse the test tubes from part A and bring the hot water bath to a boil.
2. Add about 2 mL of appropriate sugar solutions to the respective labeled test tubes, including the water blank as described in A-1 above.

**Caution: Barfoed's reagent is caustic, rinse thoroughly with water if you get this solution on your skin or clothing.**

3. Add 2 mL of Barfoed's reagent to each tube and mix well.
4. Place all the tubes in the boiling water bath at the same time and heat for 2 minutes after the water begins to boil again.
5. Record observations on the report sheet. Do you see a red precipitate at the bottom of the tubes?
6. What conclusions can you make about your unknown now?

Discard these solutions in the sink and rinse with plenty of water.

**Part C: Test for ketohexoses with Seliwanoff's reagent.**

1. Empty, wash and rinse the test tubes from part B and make sure you have sufficient water in your boiling water bath.

**Caution: Seliwanoff's reagent is caustic, rinse thoroughly with water if you get this solution on your skin or clothing.**

2. Add about 3 mL of Seliwanoff's reagent to each labeled test tube.
3. Add 1 drop of the respective sugar solution and 1 drop of water blank to the appropriate test tubes as described in part A-1 above and mix well.

4. Place all the test tubes in the boiling water bath at the same time and heat for 3 min after the water begins to boil again.
5. Record observations on the report sheet.
6. What conclusions can you make about your unknown now?

Discard these solutions in the sink and rinse with plenty of water.

#### **Part D: Test for pentoses with Bial's reagent.**

**CAUTION!!! Be extremely careful when heating these tubes. Bial's reagent contains concentrated hydrochloric acid. Be sure to heat the tube slowly, tilting it slightly and heating along the entire side of the tube. Do not point the tube toward yourself or any of your fellow lab workers.**

1. Empty, wash and rinse the test tubes from part C and remove the flame from the hot water bath. You will use the flame to heat the test tubes directly in this part.
2. Add about 2 mL of 1% xylose, glucose, fructose and your unknown solution to their respective labeled test tube.
3. Add 3 mL of Bial's reagent to each tube and mix well.
4. Carefully heat each tube (with some agitation) directly over the burner flame. Hold the tube at a diagonal and heat along the sides of the tube rather than at the bottom to prevent eruption of the liquid from the tube. Move the tube diagonally in and out of the flame, until the mixture just begins to boil. Stop heating when the mixture begins to boil.
5. Record observations on the report sheet.
6. What conclusions can you make about your unknown now?

Discard these solutions in the sink and rinse with plenty of water.

#### **Part E: Test for starches with iodine.**

1. Empty, wash and rinse the test tubes from part D. Place the flame under the water bath; it will be needed for part F.
2. Add about 2 mL of each carbohydrate test solution or deionized water to the respective labeled test tube as described in part A-1.
3. Add 1 drop of iodine reagent to each tube and mix well.

4. Record observations on the report sheet.
5. Summarize your conclusion about the identity of your unknown carbohydrate indicating the support for this conclusion as a result of each of the tests you performed.

Discard these solutions in the sink and rinse with plenty of water.

**Part F: Enzymatic Digestion of Starch.**

1. Collect about 2 mL of saliva in a clean test tube.
2. Add 6 mL of neutral (pH 7) buffer solution to the saliva and mix well.
3. Pour half of the saliva mixture into another clean test tube and label these test tubes S1 and S2 (for saliva)
4. Add about 4 mL of neutral buffer (no saliva) to each of 2 other clean test tubes and label these B1 and B2 (for buffer).
5. Add about 10 drops of starch solution to each of the 4 tubes (2 with saliva and 2 with buffer only) and mix well. Allow the mixtures to stand while you proceed with parts A thru E.
6. After allowing the solutions to stand for at least 30 min, test for glucose with Benedict's reagent and starch with iodine as described below.
7. Add 3 drops of Benedict's reagent to one tube containing saliva (S1) and to one tube containing buffer (B1), mix well and place both tubes in the hot water bath for several minutes.
8. Add 1 drop of iodine solution to the other tube containing saliva (S2) and the other tube containing buffer (B2) and mix well. It is not necessary to heat these 2 tubes.
9. Record your observations for these 4 tubes on the report sheet and answer the questions for this part.





Name \_\_\_\_\_

Section \_\_\_\_\_

## The Chemistry of Carbohydrates

Experiment #7

Data &amp; Report Sheet

Observations for Parts A thru E.

Unknown number \_\_\_\_\_

Carbo- hydrate	Benedict Test	Barfoed Test	Seliwanoff Test	Bial Test	Iodine Test
Glucose					
Fructose					
Xylose					
Lactose					
Sucrose					
Starch					
Unknown #					
Water					

A-1. What conclusions can you make regarding your unknown after the test with Benedict's reagent? Is your unknown a reducing sugar?

B-1. What conclusions can you make about your unknown after the test with Barfoed's reagent? Is your unknown a monosaccharide? An aldose?

C-1. What conclusions can you make about your unknown after the test with Seliwanoff's reagent? How does the final color of your unknown in the Seliwanoff test support your conclusion?

D-1. What conclusions can you draw about your unknown after the test with Bial's reagent? How does the color of your unknown in the Bial test support your conclusion?

E-1. What conclusions can you make about your unknown after the test with the Iodine reagent?

Summarize your conclusions. What is the identity of your unknown (aldose or ketose; pentose or hexose; monosaccharide, disaccharide or starch)?

Name \_\_\_\_\_

Section \_\_\_\_\_

Part F. Observations. Indicate whether each test was positive or negative and indicate what color each solution was after performing the tests.

Starch mixture	Benedict's Test	Iodine Test
Starch/Solution	S1	S2
Starch/Buffer	B1	B2

F-1. Describe what caused the observed results in the Benedict's Test for starch/saliva vs. starch/buffer mixtures.

F-2. Describe what caused the observed results in the iodine test for the starch/saliva vs. starch/buffer mixtures.

F-3. What conclusions can you make regarding the action of saliva on starch?

F-4. What enzyme is involved? You may want to consult the text book for this.

F-5. What is (are) the product(s) of the reaction of starch with this enzyme?